## Installation av Philips SpeechMagic 6.1 (endast usb-drivrutiner)

1. Placera katalogen "PhilipsSpeechMagic61" på lämpligt lokalt ställe.

2. Kör filen "PhilipsSpeechMagic61\setup.exe". Bilden nedan kommer att visas. Välj Next.



3. Bilden nedan kommer att visas. Markera "I accept..." och välj Next.

License Agreement	InstallSineta wizaru	
Please read the following license	e agreement carefully.	
SOFTWARE DEVELOPMENT KIT AG SOFTWARE	GREEMENT / PHILIPS SPEECH RECO	GNITION SYSTEMS
IMPORTANT - READ CAREFULLY : I terms and conditions. This Softwar agreement between you (either an PSRS ("Philips") regarding use of th includes computer software with ac "online" or electronic documentatio installing, copying or otherwise usir this Agreement, you agree to be b	Before installing this software care re Development Kit Agreement ("Ag n individual or a single entity) and P he Philips software accompanying t iccompanying data ("SOFTWARE") ( on (together "SOFTWARE PRODUCT) ing the SOFTWARE PRODUCT, by a bound by these terms. If you do no	fully read the following reement") is a binding 'hilips Austria GmbH, his Agreement, which and may also include T"). In downloading, accepting the terms of t agree with the terms of
I accept the terms in the license	e agreement	
] I do not accept the terms in the	license agreement	
a det a la		
stalbhield		

4. Bilden nedan kommer att visas. Bocka för drivrutiner enligt nedan och välj Next.

🙀 SpeechMagic Recorder - InstallShield Wizard	
<b>Custom Setup</b> Select the program features you want installed.	
Click on an icon in the list below to change how a feature is in Sound Recorder application AudioWizard SmRecord ActiveX Control SmAudio COM Control Control Device Drivers SpeechMike Pro or SpeechMike SpeechMike Classic Foot Control with serial port adapter	stalled. Feature Description This feature installs the driver for the Philips Foot Control LFH 0210 with serial port adapter. This feature requires 86KB on your hard drive.
Install to: C:\Program Files\SpeechMagic\Recorder\6.1.0\ InstallShield Help Space < Back	Change Next > Cancel

5. Bilden nedan kommer att visas. Välj Install.



6. Bilden nedan kommer att visas. Välj Finish.



< Back	Finish	Cancel

7. Bilden nedan kommer att visas. Om du vill starta om datorn, välj Ja.



## Installation av SpeechMagic6.1 servicepack

- 1. Placera katalogen "SpeechMagic6.1 servicepack" på lämpligt lokalt ställe.
- 2. Kör filen "...\SpeechMagic6.1 servicepack\Recorder61B543SP02\setup.exe". Bilden nedan kommer att visas. Välj Yes.

Speechl	Magic Recorder
2	This setup will perform an upgrade of 'SpeechMagic Recorder'. Do you want to continue?
	Yes No

3. Följande bild kommer att visas. Välj Next.



4. Följande bild kommer att visas. Vänta...

Speech/	Magic Recorder - InstallShield Wizard
17	Please wait while the InstallShield Wizard installs SpeechMagic Recorder. This may take several minutes.
	Status:
	Validating install
1.00	
InstallShield	< Back Next > Cancel

5. Till slut visas följande bild. Välj Finish.



6. Bilden nedan kommer att visas. Om du vill starta om datorn, välj Ja.

